

Summary of laboratory means

sample 1

Laboratory	n-Hexane	Z score	n-Octane	Z score	n-Propyl acetate	Z score
Unit	mg/tube		mg/tube		mg/tube	
–	–	--	–	--	–	--
38	0,681	0,313	0,413	0,650	0,210	0,384
68	0,656	-0,066	0,398	0,263	0,196	-0,308
79	0,661	0,004	0,388	-0,005		
82	0,689	0,434	0,395	0,186	0,218	0,779
85	0,604	-0,853				
118	0,630	-0,466	0,370	-0,469	0,200	-0,101
125						
128						
131	0,625	-0,535	0,362	-0,665	0,207	0,236
135	0,610	-0,764	0,358	-0,756	0,197	-0,244
162	0,674	0,208	0,403	0,402	0,179	-1,144
195	0,808	2,230 E	0,375	-0,327	0,089	-5,584 BE
202	0,661	0,013	0,391	0,088	0,240	1,867
224	0,681	0,313	0,413	0,650	0,206	0,186
262	0,655	-0,081	0,393	0,134	0,195	-0,358
272	0,611	-0,747	0,382	-0,150	0,176	-1,297
–	–	--	–	--	–	--
Method	ISO 5725		ISO 5725		ISO 5725	
Assessment	Z <=2,000		Z <=2,000		Z <=2,000	
Mean	0,660		0,388		0,202	
Reproducibility s.d.	0,051		0,018		0,018	
Rel. reproducibility s.d.	7,74 %		4,56 %		8,73 %	
Reference value	0,658		0,397		0,214	
Target s.d.	0,066		0,039		0,020	

Laboratory	n-Hexane Z score	n-Octane Z score	n-Propyl acetate Z score
Rel. target s.d.	10,00 %	10,00 %	10,00 %
Lower limit of tolerance	0,528	0,310	0,162
Upper limit of tolerance	0,792	0,465	0,243
Type B outliers	0	0	1
Type E outliers	1	0	1
Type F outliers	0	0	0
No. of laboratories after elimination of outliers type A-D and F (without laboratories that only gave states but no measured values)	14	13	11
Explanation of outlier types			
A: Single outlier			
B: Differing laboratory mean			
C: excessive laboratory s.d.			
D: Excluded manually			
E: score outside tolerance limits			
F: Score >3,5			

Summary of laboratory means

sample 2

Laboratory	o-Xylene	Z score	n-Heptane	Z score	n-Octane	Z score	Ethylbenzene	Z score
Unit	mg/tube		mg/tube		mg/tube		mg/tube	
–	–	--	–	--	–	--	–	--
38	0,080	0,214	0,363	-0,384	0,393	0,139	0,139	0,362
68	0,080	0,214	0,387	0,251	0,392	0,113	0,142	0,586
79	0,087	1,095	0,391	0,344	0,391	0,085	0,137	0,191
82	0,076	-0,297	0,390	0,331	0,392	0,113	0,119	-1,129
85	0,084	0,725					0,141	0,511
118	0,073	-0,731	0,358	-0,511	0,348	-1,014	0,124	-0,734
125	0,086	0,980					0,144	0,735
128	0,071	-0,935					0,129	-0,383
131	0,077	-0,169	0,365	-0,331	0,356	-0,816	0,127	-0,533
135	0,075	-0,466	0,351	-0,710	0,354	-0,860	0,128	-0,488
162	0,073	-0,680	0,388	0,270	0,392	0,103	0,130	-0,294
195	0,068	-1,292	0,381	0,098	0,397	0,247	0,131	-0,242
202	0,085	0,827	0,402	0,657	0,413	0,650	0,142	0,556
224	0,079	0,086	0,386	0,225	0,383	-0,119	0,134	-0,011
262	0,083	0,597	0,364	-0,358	0,433	1,171	0,140	0,437
272	0,077	-0,169	0,382	0,119	0,395	0,190	0,140	0,437
–	–	--	–	--	–	--	–	--
Method	ISO 5725		ISO 5725		ISO 5725		ISO 5725	
Assessment	Z <=2,000		Z <=2,000		Z <=2,000		Z <=2,000	
Mean	0,078		0,378		0,388		0,134	
Reproducibility s.d.	0,006		0,015		0,023		0,007	
Rel. reproducibility s.d.	7,18 %		4,10 %		6,06 %		5,55 %	
Reference value	0,082		0,379		0,383		0,139	
Target s.d.	0,008		0,038		0,039		0,013	

Laboratory	o-Xylene Z score	n-Heptane Z score	n-Octane Z score	Ethylbenzene Z score
Rel. target s.d.	10,00 %	10,00 %	10,00 %	10,00 %
Lower limit of tolerance	0,063	0,302	0,310	0,107
Upper limit of tolerance	0,094	0,453	0,465	0,161
Type B outliers	0	0	0	0
Type E outliers	0	0	0	0
Type F outliers	0	0	0	0
No. of laboratories after elimination of outliers type A-D and F (without laboratories that only gave states but no measured values)	16	13	13	16
Explanation of outlier types				
A: Single outlier				
B: Differing laboratory mean				
C: excessive laboratory s.d.				
D: Excluded manually				
E: score outside tolerance limits				
F: Score >3,5				

Summary of laboratory means


sample 3

Laboratory	o-Xylene	Z score	Toluene	Z score	Ethylbenzene	Z score
Unit	mg/tube		mg/tube		mg/tube	
–	–	--	–	--	–	--
38	0,140	0,746	0,696	0,209	0,144	0,635
68	0,133	0,208	0,685	0,048	0,143	0,561
79	0,139	0,661	0,701	0,285	0,136	0,073
82	0,125	-0,406	0,670	-0,172	0,130	-0,399
85	0,139	0,669	0,697	0,224	0,144	0,635
118	0,126	-0,329	0,663	-0,278	0,129	-0,503
125	0,148	1,360	0,698	0,238	0,150	1,078
128	0,114	-1,250	0,666	-0,231	0,133	-0,178
131	0,128	-0,176	0,663	-0,275	0,128	-0,547
135	0,125	-0,421	0,651	-0,454	0,129	-0,488
162	0,119	-0,828	0,665	-0,250	0,128	-0,532
195	0,119	-0,843	0,726	0,652	0,130	-0,370
202	0,137	0,500	0,693	0,160	0,139	0,243
224	0,135	0,362	0,680	-0,026	0,138	0,192
262	0,127	-0,252	0,673	-0,128	0,130	-0,399
272	0,044	-6,623 BE	0,315	-5,380 BE	0,048	-6,455 BE
–	–	--	–	--	–	--
Method	ISO 5725		ISO 5725		ISO 5725	
Assessment	Z <=2,000		Z <=2,000		Z <=2,000	
Mean	0,130		0,682		0,135	
Reproducibility s.d.	0,009		0,020		0,007	
Rel. reproducibility s.d.	7,18 %		2,94 %		5,31 %	
Reference value	0,134		0,706		0,138	
Target s.d.	0,013		0,068		0,014	

Laboratory	o-Xylene Z score	Toluene Z score	Ethylbenzene Z score
Rel. target s.d.	10,00 %	10,00 %	10,00 %
Lower limit of tolerance	0,104	0,545	0,108
Upper limit of tolerance	0,156	0,818	0,162
Type B outliers	1	1	1
Type E outliers	1	1	1
Type F outliers	0	0	0
No. of laboratories after elimination of outliers type A-D and F (without laboratories that only gave states but no measured values)	15	15	15
Explanation of outlier types			
A: Single outlier			
B: Differing laboratory mean			
C: excessive laboratory s.d.			
D: Excluded manually			
E: score outside tolerance limits			
F: Score >3,5			


Summary results

Measurand: n-Hexane Mean: 0,660 mg/tube
Sample: sample 1 Reproducibility s.d.: 0,051 mg/tube
Method: ISO 5725 Rel. reproducibility s.d.: 7,74%
No. of laboratories: 14 Tolerance limits: 0,528 - 0,792 mg/tube ($|Z \text{ score}| < 2,00$)




Summary results

Measurand: n-Octane Mean: 0,388 mg/tube
Sample: sample 1 Reproducibility s.d.: 0,018 mg/tube
Method: ISO 5725 Rel. reproducibility s.d.: 4,56%
No. of laboratories: 13 Tolerance limits: 0,310 - 0,465 mg/tube ($|Z \text{ score}| < 2,00$)




Summary results

Measurand: n-Propyl acetate Mean: 0,202 mg/tube
Sample: sample 1 Reproducibility s.d.: 0,018 mg/tube
Method: ISO 5725 Rel. reproducibility s.d.: 8,73%
No. of laboratories: 11 Tolerance limits: 0,162 - 0,243 mg/tube ($|Z \text{ score}| < 2,00$)




Sample chart of Z scores

Sample: 1




Summary results

Measurand:	o-Xylene	Mean:	0,078 mg/tube
Sample:	sample 2	Reproducibility s.d.:	0,006 mg/tube
Method:	ISO 5725	Rel. reproducibility s.d.:	7,18%
No. of laboratories:	16	Tolerance limits:	0,063 - 0,094 mg/tube ($ Z \text{ score} < 2,00$)




Summary results

Measurand: n-Heptane Mean: 0,378 mg/tube
Sample: sample 2 Reproducibility s.d.: 0,015 mg/tube
Method: ISO 5725 Rel. reproducibility s.d.: 4,10%
No. of laboratories: 13 Tolerance limits: 0,302 - 0,453 mg/tube ($|Z \text{ score}| < 2,00$)




Summary results

Measurand: n-Octane Mean: 0,388 mg/tube
Sample: sample 2 Reproducibility s.d.: 0,023 mg/tube
Method: ISO 5725 Rel. reproducibility s.d.: 6,06%
No. of laboratories: 13 Tolerance limits: 0,310 - 0,465 mg/tube ($|Z \text{ score}| < 2,00$)




Summary results

Measurand:	Ethylbenzene	Mean:	0,134 mg/tube
Sample:	sample 2	Reproducibility s.d.:	0,007 mg/tube
Method:	ISO 5725	Rel. reproducibility s.d.:	5,55%
No. of laboratories:	16	Tolerance limits:	0,107 - 0,161 mg/tube ($ Z \text{ score} < 2,00$)




Sample chart of Z scores

Sample: 2




Summary results

Measurand:	o-Xylene	Mean:	0,130 mg/tube
Sample:	sample 3	Reproducibility s.d.:	0,009 mg/tube
Method:	ISO 5725	Rel. reproducibility s.d.:	7,18%
No. of laboratories:	15	Tolerance limits:	0,104 - 0,156 mg/tube ($ Z \text{ score} < 2,00$)




Summary results

Measurand:	Toluene	Mean:	0,682 mg/tube
Sample:	sample 3	Reproducibility s.d.:	0,020 mg/tube
Method:	ISO 5725	Rel. reproducibility s.d.:	2,94%
No. of laboratories:	15	Tolerance limits:	0,545 - 0,818 mg/tube ($ Z \text{ score} < 2,00$)




Summary results

Measurand: Ethylbenzene Mean: 0,135 mg/tube
Sample: sample 3 Reproducibility s.d.: 0,007 mg/tube
Method: ISO 5725 Rel. reproducibility s.d.: 5,31%
No. of laboratories: 15 Tolerance limits: 0,108 - 0,162 mg/tube ($|Z \text{ score}| < 2,00$)



Sample chart of Z scores

Sample: 3



Questions and Answers

Participant	kind of tube (NIOSH or TYPE B)	analytical method	solvent	volume	carrier gas
38	Dräger Typ NIOSH	Hausmethode	Schwefelkohlenstoff + 2% Methanol	4 ml	Helium
68	NIOSH	GC/FID	CS2	1 ml	Helium
79	NIOSH	BIA 6265	CS2	0,5 ml	Wasserstoff
82	NIOSH	ja	CS2 / Isopropanol (80/20)	1 ml	Helium
85	NIOSH	Ja	CS2	1ml	Stickstoff
118	NIOSH	in Anlehnung an die VDI 2100 Blatt 2	CS2	2mL	Stickstoff
125	NIOSH	NIOSH Nr. 1022	Schwefelkohlenstoff	1 ml	Helium
128	Niosh	Niosh 1501: 2003	Carbon disulfide	1ml	helium
131	NIOSH	own, based on NIOSH and OSHA methods	carbon disulfide	2ml	helium
135	Niosh	Hausmethode	Schwefelkohlenstoff/Methanol , 100/1	10 ml	Helium
162	NIOSH	Hausmethode	CS2	1 ml	Wasserstoff
195	NIOSH	internal method	CARBON DISULFIDE	1 mL	He
202	NIOSH	NIOSH 1501	carbon sulfide	1mL	helium
262	NIOSH	VDI 3865 Bl.3	CS2	1 ml	Helium
272	NIOSH	DFG 3	Temäres Gemisch	5 ml	Wasserstoff

Participant	injection	autosampler	kind of injection	analytical column	detector
38	Split	Ja	Flüssiginjektion	Varian VF1-MS 60m x 0,32mm x 1µm	MSD 5975 C
68	split	Ja	Lösung	Vocol von Supelco	FID
79	split	nein	Lösung	CP Sil 5 CB	FID
82	ja	ja	split	Agilent HP-5 30m X 0,32mm X 0,25µm	FID
85	Autosampler	Ja	onColumn	DB5-DBWax	FID
118	on-column	ja	Lösung	CP Sil 5 CB / CP-Wax 57 CB	FID
125	1 µl	Agilent 6890	Split	Varian Factor 4 624 MS, 20 m x 0,15 mm ID x 0,84 µm FD	MSD
128	2 microlitri	yes	split	ZB-Wax (30 m - 0.53 mm - 1.0 µm)	FID
131	split 1:10	yes	solvent	Agilent HP-5, HP-InnoWax	GC-FID

Ringversuch organic solvents NIOSH 2012

Participant	injection	autosampler	kind of injection	analytical column	detector
135	Splitlos	Ja	Lösung	Agilent DB1 und Restek RTX 200	FID
162	Split	Ja	Lösung	CP Sil PONA CB 50 m x 0,21 mm ID x 0,5 µm Film	FID
195	1 uL	atosystem-xl perkin elmer	split	élite-db5	fid
202	split	yes	solution	DB624	FID
262	split 1:10	CTC	flüssig	Varian VF 624 ms 60 m	MS
272	1 µl	ja	splitlos	DB5	FID

Participant	data evaluation	date of analysis
38	Interne Standardisierung	15./16.03.2012
68	interner Standard	20.3.2012
79	interner Standard	13. KW 2012
82	Interner Standard	2012-03-14
85	interner Standard	12.03.2012 & 13.03.2012
118	interner Standard	22.03.2012
125	Agilent Enhanced Chemstation, interner Standard Toluol - d8	23.03.2012
128		12/03/2012
131	external standard	08.03.2012
135	externer Standard	17.03.2012
162	Interner Standard	01.04.2012
195	lims	23/03/2012
202	external standard	20/03/2012
262	Flächen-Integration, MS Workstation	27.03.2012
272	externer Standard	26.03.12