

# **Round-robin tests for in-house measuring laboratories**

Institut für Arbeitsschutz der Deutschen Gesetzlichen Unfallversicherung (IFA)  
Institute for Occupational Safety and Health of the German Social Accident Insurance

B. Maybaum, K. Gusbeth, Prof. Dr. D. Breuer  
Alte Heerstraße 111, 53757 Sankt Augustin  
proficiency-testing@dguv.de, +49 2241 231 2549

## ***Results and Evaluation***

Organic solvents with sampling

03/04 November 2015

# Summary of laboratory test results

Sample 1

	4-Methyl-2-Pentanone Z score		2-Butanone Z score		Cyclopentanone Z score	
Unit	mg/m <sup>3</sup>		mg/m <sup>3</sup>		mg/m <sup>3</sup>	
5	64,56	0,65	134,09	0,50	36,08	-1,97
18	62,10	0,24	125,50	-0,17	48,70	0,84
37	66,30	0,94	128,40	0,05	25,30	-4,37 FE
40	65,05	0,73	127,40	-0,03	47,84	0,65
72	67,70	1,17	143,40	1,23	23,60	-4,74 FE
78	58,45	-0,36	118,79	-0,70	46,29	0,31
111	58,20	-0,40	119,10	-0,68	46,20	0,29
146	54,00	-1,09	125,00	-0,21	23,30	-4,81 FE
158	59,80	-0,14	135,30	0,59	48,70	0,84
236	60,04	-0,10	131,98	0,33	45,56	0,14
271	50,70	-1,64	116,00	-0,92	39,90	-1,12
281			160,00	2,53 BE		
-	-	--	-	--	-	--
Method	ISO 5725-2		ISO 5725-2		ISO 5725-2	
Assessment	Z <=2,00		Z <=2,00		Z <=2,00	
No. of laboratories that submitted results	11		12		11	
Mean	60,63		127,72		44,91	
Reproducibility s.d.	5,23		8,16		4,54	
Rel. reproducibility s.d.	8,63 %		6,39 %		10,11 %	
Reference value	60,90		122,90		47,20	
Target s.d.	6,06		12,77		4,49	
Rel. SDPA	10,00 %		10,00 %		10,00 %	
Lower limit of tolerance	48,50		102,18		35,93	
Upper limit of tolerance	72,75		153,27		53,89	
Type B outliers			1			
Type F outliers					3	
No. of laboratories after elimination of outliers type A-D and F (without laboratories that only gave states but no measured values)	11		11		8	
Explanation of outlier types						
A: Single outlier	Grubbs					
B: Differing laboratory mean	Grubbs					
C: Excessive laboratory s.d.	Cochran					
D: Excluded manually						
E: mean outside tolerance limits						
F:  Z-Score >3,5						
L: Differing laboratory mean (Grubbs II)	Grubbs II					

## Summary of laboratory test results

Sample 2

	n-Butyl acetate Z score		Toluene Z score		m-Xylene Z score		n-Hexane Z score		n-Octane Z score		Ethylbenzene Z score	
Unit	mg/m <sup>3</sup>		mg/m <sup>3</sup>		mg/m <sup>3</sup>		mg/m <sup>3</sup>		mg/m <sup>3</sup>		mg/m <sup>3</sup>	
5	76,63	0,51	64,96	0,42	60,40	-4,26 BE	58,85	0,90	131,76	0,42	94,10	8,62 BE
18	80,30	1,01	59,20	-0,50	100,90	-0,41	52,60	-0,26	134,00	0,59	50,10	-0,09
37	78,30	0,74	69,70	1,18	112,50	0,70	71,90	3,32 E	152,10	2,02 E	59,80	1,83
40	77,29	0,60	65,42	0,50	110,71	0,53	49,66	-0,80	131,36	0,38	52,63	0,41
72	71,90	-0,14	64,30	0,32	103,30	-0,18	53,70	-0,05	133,40	0,55	50,40	-0,03
78	73,56	0,09	61,36	-0,15	106,83	0,16	48,35	-1,04	127,08	0,05	50,82	0,06
111	66,60	-0,87	56,60	-0,92	96,20	-0,85	43,00	-2,03 E	113,60	-1,02	45,10	-1,08
146	73,00	0,01	63,30	0,16	106,00	0,08	47,30	-1,24	127,00	0,04	49,80	-0,15
158	60,30	-1,73	57,80	-0,72	97,50	-0,73	53,10	-0,16	105,70	-1,64	45,70	-0,96
236	66,32	-0,91	56,46	-0,94	111,81	0,63	60,26	1,16	117,53	-0,71	45,59	-0,98
271	68,00	-0,68	55,70	-1,06	89,10	-1,53	51,00	-0,55	118,00	-0,67	46,00	-0,90
281	83,00	1,38	73,00	1,71	122,00	1,60	58,00	0,75	179,00	4,15 BE	60,00	1,87
–	–	--	–	--	–	--	–	--	–	--	–	--
Method	ISO 5725-2		ISO 5725-2		ISO 5725-2		ISO 5725-2		ISO 5725-2		ISO 5725-2	
Assessment	Z <=2,00		Z <=2,00		Z <=2,00		Z <=2,00		Z <=2,00		Z <=2,00	
No. of laboratories that submitted results	12		12		12		12		12		12	
Mean	72,93		62,32		105,17		53,98		126,50		50,54	
Reproducibility s.d.	6,65		5,49		9,13		7,55		12,50		5,28	
Rel. reproducibility s.d.	9,12 %		8,82 %		8,69 %		13,99 %		9,88 %		10,44 %	
Reference value	70,50		61,30		90,60		50,20		128,80		45,30	
Target s.d.	7,29		6,23		10,52		5,40		12,65		5,05	
Rel. SDPA	10,00 %		10,00 %		10,00 %		10,00 %		10,00 %		10,00 %	
Lower limit of tolerance	58,35		49,85		84,13		43,18		101,20		40,43	
Upper limit of tolerance	87,52		74,78		126,20		64,77		151,80		60,65	
Type B outliers					1				1		1	
No. of laboratories after elimination of outliers type A-D and F (without laboratories that only gave states	12		12		11		12		11		11	

	n-Butyl acetate Z score	Toluene Z score	m-Xylene Z score	n-Hexane Z score	n-Octane Z score	Ethylbenzene Z score
but no measured values)						
Explanation of outlier types						
A: Single outlier	Grubbs					
B: Differing laboratory mean	Grubbs					
C: Excessive laboratory s.d.	Cochran					
D: Excluded manually						
E: mean outside tolerance limits						
F: $ Z\text{-Score}  > 3,5$						
L: Differing laboratory mean (Grubbs II)	Grubbs II					
	1-Methoxy-2-propanole Z score	1-methoxy-2-propyl acetate Z score				
Unit	mg/m <sup>3</sup>	mg/m <sup>3</sup>				
5	137,03	1,91	114,04	5,87	BE	
18	109,60	-0,47	77,20	0,75		
37	59,90	-4,79	73,90	0,29	BE	
40	121,10	0,52	78,04	0,86		
72	123,60	0,74	87,20	2,14	E	
78	114,30	-0,07	73,32	0,21		
111	114,10	-0,08	64,40	-1,04		
146	102,00	-1,14	73,90	0,29		
158	126,20	0,97	62,50	-1,30		
236	107,72	-0,64	56,90	-2,08	E	
271	95,00	-1,74	71,00	-0,12		
-	-	--	-	--		
Method	ISO 5725-2		ISO 5725-2			
Assessment	$ Z  \leq 2,00$		$ Z  \leq 2,00$			
No. of laboratories that submitted results	11		11			
Mean	115,06		71,84			
Reproducibility s.d.	12,36		8,69			

---

	1-Methoxy-2-propanole Z score	1-methoxy-2-propyl acetate Z score
Rel. reproducibility s.d.	10,74 %	12,10 %
Reference value	104,60	69,30
Target s.d.	11,51	7,18
Rel. SDPA	10,00 %	10,00 %
Lower limit of tolerance	92,05	57,47
Upper limit of tolerance	138,08	86,20
Type B outliers	1	1
No. of laboratories after elimination of outliers type A-D and F (without laboratories that only gave states but no measured values)	10	10

---

# Summary of laboratory test results

Sample 3

## Benzene Z score

Unit	$\mu\text{g}/\text{m}^3$	
5	140,3000	1,16
18	121,1000	-0,36
37	< 1000,0000	
40	124,0500	-0,13
72	130,0000	0,34
78	132,8900	0,57
111	128,0000	0,18
146	117,0000	-0,69
158	118,4000	-0,58
236	111,0200	-1,17
271	134,0000	0,66


Method	ISO 5725-2
Assessment	$ Z  \leq 2,00$
No. of laboratories that submitted results	11
Mean	125,6760
Reproducibility s.d.	8,9810
Rel. reproducibility s.d.	7,15 %
Reference value	127,4000
Target s.d.	12,5676
Rel. SDPA	10,00 %
Lower limit of tolerance	100,5408
Upper limit of tolerance	150,8112
No. of laboratories after elimination of outliers type A-D and F (without laboratories that only gave states but no measured values)	10

### Explanation of outlier types

A: Single outlier	Grubbs
B: Differing laboratory mean	Grubbs
C: Excessive laboratory s.d.	Cochran
D: Excluded manually	
E: mean outside tolerance limits	
F: $ Z\text{-Score}  > 3,5$	
L: Differing laboratory mean (Grubbs II)	Grubbs II


## Summary results

Measurand:	4-Methyl-2-Pentanone	Mean:	60,63 mg/m <sup>3</sup>
Sample:	1	Reproducibility s.d.:	5,23 mg/m <sup>3</sup>
Method:	ISO 5725-2	Relative reproducibility s.d.:	8,63%
No. of laboratories:	11	Reference value:	60,90 mg/m <sup>3</sup>
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	48,50 - 72,75 mg/m <sup>3</sup> ( Z-Score  <= 2,00)



## Summary results


Measurand:	2-Butanone	Mean:	127,72 mg/m <sup>3</sup>
Sample:	1	Reproducibility s.d.:	8,16 mg/m <sup>3</sup>
Method:	ISO 5725-2	Relative reproducibility s.d.:	6,39%
No. of laboratories:	11	Reference value:	122,90 mg/m <sup>3</sup>
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	102,18 - 153,27 mg/m <sup>3</sup> ( Z-Score  <= 2,00)






## Summary results

Measurand:	Cyclopentanone	Mean:	44,91 mg/m <sup>3</sup>
Sample:	1	Reproducibility s.d.:	4,54 mg/m <sup>3</sup>
Method:	ISO 5725-2	Relative reproducibility s.d.:	10,11%
No. of laboratories:	8	Reference value:	47,20 mg/m <sup>3</sup>
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	35,93 - 53,89 mg/m <sup>3</sup> ( Z-Score  <= 2,00)




# Sample chart of Z-scores

Sample 1




## Summary results

Measurand:	n-Butyl acetate	Mean:	72,93 mg/m <sup>3</sup>
Sample:	2	Reproducibility s.d.:	6,65 mg/m <sup>3</sup>
Method:	ISO 5725-2	Relative reproducibility s.d.:	9,12%
No. of laboratories:	12	Reference value:	70,50 mg/m <sup>3</sup>
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	58,35 - 87,52 mg/m <sup>3</sup> ( Z-Score  <= 2,00)




## Summary results

Measurand:	Toluene	Mean:	62,32 mg/m <sup>3</sup>
Sample:	2	Reproducibility s.d.:	5,49 mg/m <sup>3</sup>
Method:	ISO 5725-2	Relative reproducibility s.d.:	8,82%
No. of laboratories:	12	Reference value:	61,30 mg/m <sup>3</sup>
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	49,85 - 74,78 mg/m <sup>3</sup> ( Z-Score  <= 2,00)




## Summary results

Measurand:	m-Xylene	Mean:	105,17 mg/m <sup>3</sup>
Sample:	2	Reproducibility s.d.:	9,13 mg/m <sup>3</sup>
Method:	ISO 5725-2	Relative reproducibility s.d.:	8,69%
No. of laboratories:	11	Reference value:	90,60 mg/m <sup>3</sup>
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	84,13 - 126,20 mg/m <sup>3</sup> ( Z-Score  <= 2,00)




## Summary results

Measurand:	n-Hexane	Mean:	53,98 mg/m <sup>3</sup>
Sample:	2	Reproducibility s.d.:	7,55 mg/m <sup>3</sup>
Method:	ISO 5725-2	Relative reproducibility s.d.:	13,99%
No. of laboratories:	12	Reference value:	50,20 mg/m <sup>3</sup>
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	43,18 - 64,77 mg/m <sup>3</sup> ( Z-Score  <= 2,00)




## Summary results

Measurand:	n-Octane	Mean:	126,50 mg/m <sup>3</sup>
Sample:	2	Reproducibility s.d.:	12,50 mg/m <sup>3</sup>
Method:	ISO 5725-2	Relative reproducibility s.d.:	9,88%
No. of laboratories:	11	Reference value:	128,80 mg/m <sup>3</sup>
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	101,20 - 151,80 mg/m <sup>3</sup> ( Z-Score  <= 2,00)



## Summary results


Measurand:	Ethylbenzene	Mean:	50,54 mg/m <sup>3</sup>
Sample:	2	Reproducibility s.d.:	5,28 mg/m <sup>3</sup>
Method:	ISO 5725-2	Relative reproducibility s.d.:	10,44%
No. of laboratories:	11	Reference value:	45,30 mg/m <sup>3</sup>
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	40,43 - 60,65 mg/m <sup>3</sup> ( Z-Score  <= 2,00)






## Summary results

Measurand:	1-Methoxy-2-propanol	Mean:	115,06 mg/m <sup>3</sup>
Sample:	2	Reproducibility s.d.:	12,36 mg/m <sup>3</sup>
Method:	ISO 5725-2	Relative reproducibility s.d.:	10,74%
No. of laboratories:	10	Reference value:	104,60 mg/m <sup>3</sup>
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	92,05 - 138,08 mg/m <sup>3</sup> ( Z-Score  <= 2,00)



## Summary results


Measurand:	1-Methoxy-2-propyl acetate	Mean:	71,84 mg/m <sup>3</sup>
Sample:	2	Reproducibility s.d.:	8,69 mg/m <sup>3</sup>
Method:	ISO 5725-2	Relative reproducibility s.d.:	12,10%
No. of laboratories:	10	Reference value:	69,30 mg/m <sup>3</sup>
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	57,47 - 86,20 mg/m <sup>3</sup> ( Z-Score  <= 2,00)



# Sample chart of Z-scores


Sample 2

Measurand




## Summary results

Measurand:	Benzene	Mean:	125,68 µg/m³
Sample:	3	Reproducibility s.d.:	8,98 µg/m³
Method:	ISO 5725-2	Relative reproducibility s.d.:	7,15%
No. of laboratories:	10	Reference value:	127,40 µg/m³
Relative target s.d.:	10,00% (Limited)	Range of tolerance:	100,54 - 150,81 µg/m³ ( Z-Score  <= 2,00)



# Sample chart of Z-scores

Sample 3



## Questions and Answers

Participant	Sample carrier	Kind of pump	Flow rate
5	Aktivkohle	Fa. SKC, PocketPump	70 und 250 ml/min
18	Silicagel Typ BIA und Aktivkohle Typ BIA	SG350 und SG4000ex, GSA Messgerätebau GmbH	ca. 0,083 l/min und ca. 0,333 l/min
37	A-Kohle Typ BIA und Silicagel Typ BIA, Fa. Dräger	Gilian PP1EX, LFS-113 DC	300 - 330 ml/min
40	Aktivkohle, Silicagel	Gilian LFS-113DC	15 bis 200 ml/min
72	Aktivkohle	SG 350 / SG 4000	200 bzw. 333 bzw. 1000 ml/min
78	Aktivkohle und Silicagel	Gilair und GSA (DEHA)	Gas1+2: 0,3 L/Min. und Gas 3 1,0 L/Min.
111	PG1=Silica, PG2+3=AK	SG-350	PG1=75ml/min, PG2+3=350ml/min
146	Aktivkohle Typ BIA, Silicagel Typ BIA (2-Butanon)	GilAir 5	333 mL/min. 83 mL/min.
158	Tenax TA	Gilian LFS-113 DC Low Flow Sampler	0,010 L/min und
236	Tenax TA	PSA LFS113, Gilian	10 mL/min
271	A-Kohle + Silicagel (Ketone)	GSA SG 350 + GSA SG 4000	0,33 l/min + 0,083 l/min (Ketone) + 1 l/min (Benzol)
281		SKC Delux	

Participant	Volume flow measurement	Sampling time
5	Fa. DEHA, Defender	2 h
18	Digitaler Durchflussmesser, Flowmeter 4140, TSI GmbH	120 min, gemäß IFA-Arbeitsmappe 6265, 7708, 7322, 7569, 7732, 7733
37	Analyt Massenflussmesser GFM 17, 0 - 500 ml/min, geeicht auf Luft	30 min und 60 min
40	DryCal DC-Lite	140 min
72	Defender 520 (50-5000 ml/min)	1 h bzw. 2 h
78	Gilibrator	Gas 1+2: 60 Min. und Gas 3: 120 Min.
111	Bios Defender 510M	120 min.
146	Gilibrator	90 Minuten
158	Bios Int. Corp. Defender 510-L Rev C1	60 min und 30 min
236	Gilibrator	15, 30, 60 min
271	DryCal DC Lite	2 h

Participant	Analytical method	Sample preparation
5	Hausmethode	Desorption mit CS2

### Round-robin test Organic solvents with sampling 1/2015

Participant	Analytical method	Sample preparation
18	GCMS / IFA Arbeitsmappe 6265, 7708, 7322, 7569, 7732, 7733, 7735	-
37	Hausmethode	AK + 3ml Benzylalkohol, Headspace, Silicagel + 5ml Aceton, 30min Ultraschall, Direktinj.
40	validierte Hausmethoden	Schwefelkohlenstoff 1,5 ml; Dichlormethan/Methanol (9:1) 5 ml; Dichlormethan/Methanol/Wasser (65/35/2) 5 ml
72	validierte eigene SOP in Anlehnung an IFA Arbeitsmappe	Benzylalkohol und Gemisch Dichlormethan - Schwefelkohlenstoff - Methanol
78	IFA-Arbeitsmappe	Gas 1+2: 20 mL; Gas 3: 5 mL
111	PG1 = IFA7708, PG2 = IFA7322,7732, PG3 = IFA6265	PG1+PG2 = 2ml ternäres Gemisch, PG3 = 2ml CS2
146	2-Butanon DFG Methode Nr. 6 bzw. IFA Methode 7708	
158	VOC in Luft mittels Thermodesorption	Thermodesorption
236	Hausmethode	Thermodesorption
271	IFA 7733, IFA 6265, IFA 7569, IFA 7322, IFA 7732, IFA 7708	terneres Gemisch (65% Dichlormethan, 30% Schwefelkohlenstoff und 5 % Methanol) bzw. 100% Schwefelkohlenstoff

Participant	Carrier gas	Injection	Analytical column
5	Helium	split	HP1
18	-	-	-
37	Helium 5.0	Split	60m Rtx Volatile, 60m Stabilwax DA
40	Helium	split	DBWAX; DB1; HP5-MS; RTx-35
72	Helium 6.0	split	Phenomenex Zebron ZB-WAX, 30 m x 0,25 mm x0,25 mm
78	Stickstoff	split	DB-1
111	Helium	split	SE54
158	Helium	Inletsplit und Outletsplit	Restek RTX-1 60m, ID 0,25 x 1 µm
236	He	split/split	HP-5 50m ID=0.32mm Film=1.05µm
271	Stickstoff bzw. Helium (Ketone)	Split	Macherey Nagel Optima 1; 30 m * 0,32 mm ID; 5 µm Filmdicke bzw. Thermo TG-1301 MS; 30 m * 0,25 mm ID; 0,25 µm Filmdicke (Ketone)

Participant	Detector	Data evaluation	Recovery rate	Date of analysis
5	FID	ISTD	Ja	05.11.2015
18	-	-	-	05.11.-25.11.2015
37	FID	externer Standard	nein	19.11.2015
40	FID; massenselektiver Detektor	interner Standard	Ja	05.11.2015;
72	FID	externer Standard, 4 und 6 Pkt. Kalibrierung	ja (0,7 bis 0,95)	9.11. - 17.11.2015
78	FID	interner Standard	Ja	47. und 48. KW
111	PG1 = FID, PG2+PG3 = MS	interner Standard	nein	10.-27.11.15

Round-robin test Organic solvents with sampling 1/2015

Participant	Detector	Data evaluation	Recovery rate	Date of analysis
146				25.11.2015
158	MSD Agilent 5975C	Interner Standard	Es wurde ein Kontrollstandard eingesetzt	09.11.2015 - 20.11.2015
236	FID	externe Kalibration	nein	11. - 18.11.2015
271	FID, Range = 1; At = -4 bzw. MS (Ketone)	Interner Standard (Schwefelkohlenstoff mit Chlorbenzol, tern. Gemisch mit Undecan)	nein	10.11.2015 - 24.11.2015